

High Quality Pipette Tips Now Available from IVUS



IVUS high quality pipette tips to fit your finest work

- ✓ Available with or without an aerosol barrier filter
- ✓ Available in standard and low-retention
- ✓ Superior quality to other standard and low-retention tips
- ✓ Sterile, bulk and reload options available

User-friendly rack

- ✓ Front and back snap card tab placement is ideal for all hand sizes
- ✓ Locking lid to protect the tips from bumps and spills
- ✓ Built with stability for ease of use with multi-channel pipetting
- ✓ Back hinge makes it easy to open and close with one hand
- ✓ Top can be removed for easy access
- ✓ Easy to read tip size

Reload Options



Eco-friendly attributes

- ✓ Reload system for low-retention non-filter tips
- ✓ Fully recyclable trays and racks

Se alle spidser her



Find The IVUS Tip That Meets Your Research Needs

IVUS quality tips are available in the following categories:

Racked Non Filter - Non Low Retention Non Sterile		
IV-7017	200µl Racked No Filter No Low retention	10 Racks of 96 per Pack
IV-7001	1000µl Racked No Filter No Low retention	8 Racks of 96 per Pack
Racked Non Filter - Non Low Retention Pre-Sterile		
IV-7002	1000µl Racked Non Low Retention Pre-Sterile	8 Racks of 96 per Pack
Racked Non Filter - Low Retention Pre-Sterile		
IV-7003	10xl Racked LR Pre-Sterile	10 Racks of 96 per Pack
IV-7004	200µl Racked LR Pre-Sterile	10 Racks of 96 per Pack
Racked Filter - Non Low Retention Pre-Sterile		
IV-7005	1000µl Filter Racked, Pre-Sterile	8 Racks of 96 per Pack
Racked Filter - Low Retention Pre-Sterile		
IV-7006	10µl Filter Racked LR Pre-Sterile	10 Racks of 96 per Pack
IV-7007	10xl Filter Racked LR Pre-Sterile	10 Racks of 96 per Pack
IV-7008	10E Filter Racked LR Pre-Sterile	10 Racks of 96 per Pack
IV-7009	20µl Filter Racked LR Pre-Sterile	10 Racks of 96 per Pack
IV-7010	100µl Filter Racked LR Pre-Sterile	10 Racks of 96 per Pack
IV-7011	200µl Filter Racked LR Pre-Sterile	10 Racks of 96 per Pack
IV-7012	1000µl Filter Racked LR Pre-Sterile	8 Racks of 96 per Pack
Reloads Low retention - Non Sterile		
IV-7013	200µl Reload Low Retention	10 Racks of 96 per Pack
IV-7016	10µl Reload Low Retention	10 Racks of 96 per Pack
IV-7018	1000µl Reload Low Retention	8 Racks of 96 per Pack
IV-7019	1250µl Reload Low Retention	8 Racks of 96 per Pack
IV-7020	Empty Trays for Reloads 10 and 200ul Tips	
IV-7021	Empty Trays for Reloads 1000 and 1250ul Tips	
Bulk Non Filter Tips Low Retention		
IV-7014	200ul GRAD TIP BULK LOW	1000 Tips/Pack, 4 Packs/Case
Racked Non Filter, Non Low Retention Non Sterile		
IV-7015	10µl Racked No Filter No Low retention	10 Racks of 96 per Pack



Original IVUS Quality